RIDASCREEN[®] Melengestrolacetat

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RIDASCREEN[®] Melengestrolacetat

Brief information

RIDASCREEN[®] Melengestrolacetat (Art. No.: R6502) is a competitive enzyme immunoassay for the quantitative analysis of melengestrol acetate in bovine perirenal fat and muscle meat.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit. The test kit is sufficient for 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	extraction over night centrifugation, evaporation, freezing, centrifugation chromatographic purification with RIDA [®] C18 columns (Art. No.: R2002)
Time requirement:	sample preparation (for 10 samples) approx. 18 h test implementation (incubation time) approx. 2 h 30 min
Detection limit:	bovine perirenal fat 0.3 ppb muscle meat 0.075 ppb
Recovery rate:	in spiked negative control bovine perirenal fat approx. 65 % muscle meat approx. 74 %
Cross-reactions:	The specificity of RIDASCREEN [®] Melengestrolacetat was established by analysing the cross-reactivity to corresponding substances. Melengestrol acetate (MLGA) 100.0 % Megestrol acetate (MGA) 10.0 % Medroxyprogesterone acetate (MPA) 6.6 % 6 α -Methyl-17 α -OH-Progesterone acetate 5.2 % Chlormadinone acetate (CMA) 4.8 % 17 α -Acetoxyprogesterone 3.9 % 16 α -Methylprogesterone 0.0035 % cross reaction of all other substances being checked in the assay \leq 0.003 %

1. Intended use

The RIDASCREEN® Melengestrolacetat test is a competitive enzyme immunoassay for the quantitative analysis of melengestrol acetate in bovine perirenal fat and muscle meat. The performance of the test is adjusted to the current legal demands and well suited for the control of meat exported by the USA and Canada.

2. General

Synthetic gestagens can be used for estrus inhibition or synchronisation after break off and improvement of fertility. In addition one of these substances, melengestrol acetate (MLGA), is a licensed growth promoting feed additive for heifers in the USA and Canada. The admitted dose is 0.5 mg per day and head, given as a feed premix. Its mode of anabolic action is unclear, but it stimulates the ovarian synthesis of the endogenous anabolic steroid estradiol and that may have androgenic side effects. MLGA belongs to the most active synthetic gestagens. Its oral bioactivity is about 10 or 100 times higher than the activities of the gestagens chlormadinone acetate (CMA) or medroxyprogesterone acetate (MPA), resp. If given parenterally the hormonal activity of MLGA is still 125 times higher than the one of progesterone. Due to its strongly lipophilic properties, MLGA is accumulated in fat 200-fold higher than in blood plasma.

In the European Union the use of sexual hormones for growth promoting purposes is generally forbidden since 1988, and also the import of meat from hormone treated cattle into the EU is prohibited.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-MLGA antibodies. Standards or sample solution, MLGA enzyme conjugate and anti-MLGA antibodies are added. Free and enzyme conjugated MLGA compete for the antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-MLGA antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Enzyme substrate (urea peroxide) and chromogen (tetramethylbenzidine) are added to the wells and incubated. Bound enzyme conjugate converts the colorless chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the melengestrol acetate concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with capture antibodies
- 6 x MLGA standard solutions, 1.3 ml each
0 ppb (zero standard), 0.15 ppb, 0.45 ppb, 1.35 ppb, 4.05 ppb, 12.15 ppb
MLGA in aqueous solution, ready to use
- 1 x Conjugate (0.7 ml)red cap
peroxidase conjugated MLGA
concentrate
- 1 x Anti-MLGA antibody (0.7 ml) black cap
concentrate
- 1 x Substrate (7 ml)green cap
contains urea peroxide
- 1 x Chromogen (7 ml) blue cap
contains tetramethylbenzidine
- 1 x Stop solution (14 ml)yellow cap
contains 1 N sulfuric acid
- 1 x Buffer (25 ml)
Conjugate and antibody dilution buffer

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- freezer (-60 °C / -76 °F)
- water bath (40 °C and 60 °C / 104 °F and 140 °F)
- ultra turrax
- refrigerated centrifuge
- graduated pipettes
- variable 20 µl-, 50 µl-, 100 µl-, 200 µl- and 1000 µl-micropipettes
- RIDA® C18 columns (Art. No.: R2002)

5.2. Reagents:

- petroleum ether
- methanol
- methanol/20 mM tris-HCl (20/80 (v/v)): dissolve 2.42 g tris-(hydroxymethyl) aminomethane in approx. 700 ml distilled water + 200 ml methanol (100 %), adjust pH 8.5 with approx. 5 M HCl and fill up to 1000 ml with distilled water
- 67 mM phosphate buffer, pH 7.2: 1.79 g $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$; 9.61 g $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$ + 8.7 g NaCl, fill up to 1000 ml with distilled water
- Tween 20 (Sigma, Order No.: P-2287)

6. Warnings and precautions for the users

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag and reseal them together with the desiccant provided at 2 - 8 °C (35 - 46 °F).

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Any bluish coloration of the chromogen solution solution prior to test implementation
- A value of less than 0.6 absorbance units ($A_{450 \text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

9.1. Bovine perirenal fat

- grind the frozen kidney fat by scraping with a scalpel
- weigh 1 g of sample into a 20 ml vial for extraction and add 15 ml petroleum ether and shake softly in a waterbath at 40 °C (104 °F) for extraction overnight
- centrifuge the extract : 15 min/approx. 2000 g/-15 °C (5 °F)
(if a refrigerated centrifuge is not available, freeze the sample before centrifugation)
- decant the petroleum ether supernatant into a new vial and evaporate it in a 60 °C (140 °F) waterbath
- redissolve the residue in 2 ml methanol and vortex for 20 sec. (close vial tightly)
- for degreasing: freeze 45 min at -60 °C (-76 °F)
- centrifuge: 5 min/approx. 2000 g/-15 °C (5 °F) and decant the supernatant into a new vial
- dilute the methanolic supernatant with 5 ml distilled water and purify by RIDA[®] C18 columns (Art. No.: R2002) as follows:

flow rate: 1 drop per sec

- rinse column twice with 1 ml methanol (100 %)
- equilibrate column twice with 1 ml methanol/20 mM tris-HCl, pH 8.5 (20/80 (v/v))
- apply diluted sample (approx. 7 ml)
- rinse the column twice with 1 ml methanol/20 mM tris-HCl, pH 8.5 (20/80 (v/v))
- rinse the column twice with 1 ml methanol/ distilled water (40/60 (v/v))
- remove fluid residues by positive pressure or vacuum and dry column for 2 min by floating it with air or nitrogen
- elute sample with 1 ml methanol/distilled water (80/20 (v/v)) into a new vial
- collect residues by applying positive pressure or vacuum (flow rate: 1 drop/sec)
- dilute the eluate 1 + 1 with distilled water to get a methanol/distilled water solution of 40/60 (v/v)
- employ 20 µl per well in the assay

remark:

Depending on the MLGA concentration in the sample further dilutions may be necessary. In this case use methanol/distilled water (40/60 (v/v)) for all further dilutions (i. e. always maintain 40 % methanol).

9.2. muscle meat

- cut 50 - 300 g of muscle meat into cubes and add the corresponding volume of 67 mM of phosphate buffer, pH 7.2 and homogenize with an ultra turrax
- weigh 8 g of sample into a 20 ml vial and add 8 ml petroleum ether for extraction and shake softly in a waterbath at 40 °C (104 °F) for extraction overnight
- The following sample preparation steps are identical to 9.1. !

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **MLGA enzyme conjugate** (bottle with red cap) is provided as a concentrate. Since the diluted enzyme conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the enzyme conjugate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in buffer (e. g. 200 µl conjugate concentrate + 2 ml buffer, ready to use sufficient for 4 microtiter strips).

The **anti-MLGA antibody solution** (bottle with black cap) is provided as a concentrate. Since the diluted antibody solution has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the antibody concentrate should be shaken thoroughly. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in buffer (e. g. 200 µl antibody concentrate + 2 ml buffer, ready to use sufficient for 4 microtiter strips).

For the **wash solution** use distilled water with 0.05 % tween 20 (e. g. 0.5 g tween 20 in 1000 ml of distilled water). This solution must be freshly prepared for each test series.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 20 µl of each standard solution or prepared sample to separate duplicate wells.
3. Add 50 µl of diluted enzyme conjugate to each well.
4. Add 50 µl of the diluted anti-MLGA antibody solution to each well. Mix gently by shaking the plate manually and incubate 2 h at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 250 µl of wash solution (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by rocking the plate manually and measure the absorbance at 450 nm against an air blank. Read within 30 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available to evaluate the RIDASCREEN[®] ELISA assays.

The course of the standard curve is shown in the joined Quality Assurance Certificate.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the MLGA concentration in [$\mu\text{g}/\text{kg}$].

In order to obtain the MLGA concentration in $\mu\text{g}/\text{kg}$ (ppb) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

bovine perirenal fat.....	2
muscle meat.....	0.5

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG