

Evaluation of Animal Feeds:

Carbohydrate Analysis on BEH Amide Prototypes

1/14/2010

- Sample preparation is critical for accurate quantitative analysis of carbohydrates in feed samples
- Enzymatic hydrolysis needs to be arrested for accurate determination of carbohydrates present
 - Use heat or 50% Acetonitrile extraction solvent
 - Un-arrested hydrolysis might provide useful information for total carbohydrates present
- Recommended sample preparation:
 - Extraction of fats may be necessary
 - Extract sugars with **50% ACN**
 - Homogenize for sample for 3 minutes
- Samples analyzed on BEH Amide prototypes
 - Gradient conditions necessary to elute larger polysaccharides

- First gradient method developed could NOT resolve the glucose/galactose peak pair as required by the Method Needs Statement
- Second gradient method developed utilizing linked columns
- Moderate resolution of these peaks has been obtained
 - Requires longer column (or linked columns)
 - Requires longer gradient method
 - Uses alternative mobile phase system


Carbohydrates: Modified Gradient Conditions

| Carbohydrates: Gradient Conditions | |
|---|--|
| Mobile phase | A2) 50:50 acetone/acetonitrile with 0.1% NH ₄ OH B2) 100% MilliQ Water with 0.1% NH ₄ OH |
| Gradient | 10 minute hold at 80% A2 10 minute Gradient from 80% to 40% A2 5 minute hold at 40% A2 15 minute re-equilibration in 80% A2 |
| Flow rate (mL/min) | 0.10 mL/min |
| Columns | (2 x) 2.1 x 150mm BEH Amide Prototype (1.8µm) |
| Temperature | 35°C |
| Injection | 2µL |
| Detector | Acquity ELSD |
| Gain | 200 |
| Pressure | 40 psi |
| Drift Tube | 40°C |
| Nebulizer | Cooling |
| Standard Sample | Food Sugars in 50:50 Acetonitrile/Water – 1mg/ml each (Fructose, Glucose, Sucrose, Maltose and Lactose with p-toluamide) |

Total method time: 40 minutes

Results from Previous Method

Feed Sample A




Carbohydrates: Gradient Conditions

| | |
|--------------------|---|
| Mobile phase | A2: 80% ACN with 0.2% TEA, B2: 30% ACN with 0.2% TEA |
| Gradient | 10 minute isocratic with 75% ACN, 5 minute Gradient to 50% ACN, 15 minute re-equilibration in 75% ACN |
| Flow rate (mL/min) | 0.15 mL/min |
| Column | 2.1 x 150mm BEH Amide Prototype (1.7µm) |
| Temperature | 35°C |

Results from Current Method


Feed Sample A

| Carbohydrates: Gradient Conditions | |
|------------------------------------|--|
| Mobile phase | A2: 50:50 Acetone/ACN with 0.1% NH ₄ OH, B2: 100% H ₂ O with 0.1% NH ₄ OH |
| Gradient | 10 minute isocratic with 80% A2, 10 minute Gradient to 40% A2, 5 minute hold at 40% A2, 15 minute re-equilibration in 80% A2 |
| Flow rate (mL/min) | 0.10 mL/min |
| Column | (2 x) 2.1 x 150mm BEH Amide (1.8μm) |
| Temperature | 35°C |




Increased Peak Resolution

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Spiking Experiments



Summary

- Gradient method using 0.2% TEA modified mobile phase provides excellent injection to injection and column to column reproducibility.
- Glucose/Galactose resolution is not sufficient with the this method.
- Using two 4.6 x 150mm columns with 0.1% NH₄OH modified mobile phase provides significant improvement in the resolution of the glucose/galactose pair
 - Provides better overall discretion between saccharide identities
 - Requires longer time (~ 40 minutes) than previous method
- Spiking experiments show clear distinction between glucose and galactose peaks at 1% by weight spiking level
- Using pyrrolidine as a mobile phase modifier may change selectivity and resolution enough so that a single column may be used to meet required criteria on a single 4.6 x 150mm column.

Food Sugars Analysis in Feed Samples New Data

January 18, 2010

Chromatographic Analysis of Carbohydrates

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


- Summary Highlights
 - Improved chromatographic conditions
 - Sample preparation
 - Analysis of seed Samples A-K
 - Retention time reproducibility
 - USP resolution of critical pairs

Method

ACQUITY UPLC System with ELS Detection

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Waters ACQUITY UPLC® system
- Evaporative Light Scattering Detection (ELSD)
 - Not dependent on UV absorbance
 - No derivatization is necessary
- Gradient compatible
- More sensitive than Refractive Index (RI) detection
 - At least an order of magnitude more sensitive
- Non-linear response
 - Requires a calibration curve for quantification



Detector Conditions

Waters ACQUITY ELS Detection:

| | |
|-----------------------|---|
| Gain: | 200 for sucrose and 1000 for the other 5 sugars |
| Pressure: | 40 psi |
| Drift Tube: | 40°C |
| Nebulizer: | Cooling |
| Data Rate: | 10pps |
| Filter Time Constant: | Normal |
| Data Processing: | Savitsky-Golay Smoothing (Level 17) |

Gradient Conditions Food Sugars (ELSD)

| Carbohydrates: Gradient Conditions | |
|---|--|
| Mobile phase | A) 80:20 ACN/H ₂ O with 0.1% Pyrrolidine, B) 30:70 ACN/H ₂ O with 0.1% Pyrrolidine |
| Column | 2.1 x 150 mm BEH Amide (1.7 μm) |
| Sample | Food Sugar Mixture (Fructose, Glucose, Galactose, Sucrose, Maltose and Lactose, each at 1 mg/mL in 50:50 ACN/H ₂ O) |
| Temperature | 25°C |
| Injection | 1.0μL (5μL loop with partial needle overfill) |
| Detector | Acquity ELSD |
| Detector Settings | Data rate: 10 pps, Time Constant: Normal |
| Gain | 200-300* |
| Pressure | 40 psi |
| Drift Tube | 40°C |
| Nebulizer | Cooling |
| Needle Wash | 800μL Strong needle wash of 20:80 ACN/H ₂ O, 500μL Weak needle wash of 75:25 ACN/H ₂ O |
| Seal Wash | Seal wash in 50:50 ACN/H ₂ O |
| Data Processing | Savitsky-Golay Smoothing (Level 17) |

•ELSD normalized and calibrated to Acetaminophen Test Solution (Part No. 700002387)

Gradient Conditions

The screenshot shows the Acquity Binary Solvent Manager software interface. The window title is "gr_75H5min_75_50_20min_3 in Foodsugars as System/Administrator - Instrument M...". The interface includes a menu bar (File, Edit, View, Help) and a toolbar with icons for file operations. Three instrument modules are visible: Acquity UPLC ELS Detector (ACQ-ELSD), Acquity UPLC Binary Solvent Manager (ACQ-BSM), and Acquity UPLC Sample Manager (ACQ-SM).

The main window displays the "Acquity Binary Solvent Manager" title and the "Ultra Performance LC" logo. The "General" tab is selected, showing the following settings:

- Solvents:** A2 (80:20A:H0.1%Pyrrolidir) and B2 (30:70A:H0.1%Pyrrolidir).
- Pressure Limits:** Low: 0 psi, High: 15000 psi.
- Seal Wash:** 5.0 min.

The "Gradient" section contains a table with the following data:

| | Time (min) | Flow (mL/min) | %A | %B | Curve |
|---|------------|---------------|------|------|---------|
| 1 | initial | 0.120 | 90.0 | 10.0 | initial |
| 2 | 5.00 | 0.120 | 90.0 | 10.0 | 6 |
| 3 | 25.00 | 0.120 | 40.0 | 60.0 | 6 |
| 4 | 25.50 | 0.120 | 90.0 | 10.0 | 6 |
| 5 | 45.00 | 0.120 | 90.0 | 10.0 | 6 |
| 6 | | | | | |

A "Comment:" field is located below the gradient table.

Sample Preparation


- Very simple sample preparation
- Extraction of Solid Samples
 - Weigh out sample (~3g) into 50mL centrifuge tube
 - Add 25mL of 50:50 ACN/H₂O and homogenize for at 3 minutes (mechanically)
 - Centrifuge at 4000rpm for 30 minutes
 - Collect supernatant and filter using 0.45µm PVDF syringe filter
- Depending on sample or detection method, additional sample dilutions may be necessary

Chromatograms

FEED SAMPLES (A-D) Full Scale

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


- 1. Fructose
- 2. Glucose
- 3. Galactose
- 4. Sucrose
- 5. Maltose
- 6. Lactose



FEED SAMPLES (A-D) Expanded Scale


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- 1. Fructose
- 2. Glucose
- 3. Galactose
- 4. Sucrose
- 5. Maltose
- 6. Lactose




FEED SAMPLES (E-H) Full Scale

- 1. Fructose
- 2. Glucose
- 3. Galactose
- 4. Sucrose
- 5. Maltose
- 6. Lactose



FEED SAMPLES (E-H) Expanded Scale


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



FEED SAMPLES (I-K) Full Scale


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

1. Fructose
2. Glucose
3. Galactose
4. Sucrose
5. Maltose
6. Lactose



FEED SAMPLES (I-K) Expanded Scale


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Retention Time Reproducibility

Retention Time Reproducibility (Sugars Spiked into Feed Sample K)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Each sugar was spiked at 0.1% (galactose spiked 0.15%)

Precision of Retention Times in Spiked Samples

| Name | Retention Time | %RSD |
|-----------|----------------|------|
| Fructose | 7.95 | 0.12 |
| Glucose | 9.28 | 0.18 |
| Galactose | 9.60 | 0.25 |
| Sucrose | 11.93 | 0.05 |
| Maltose | 13.28 | 0.13 |
| Lactose | 13.85 | 0.10 |

- Average retention times in the spiked samples of feed K (six replicas)
- Each sugar was spiked at 0.1% except galactose (spiked at 0.15%)

Resolution between Glucose and Galactose

| Feed Sample K | USP Resolution of Glucose & Galactose Peaks |
|----------------------|--|
| Spiked # 1 | 0.891 |
| Spiked # 2 | 0.721 |
| Spiked # 3 | 0.676 |
| Spiked # 4 | 0.829 |
| Spiked # 5 | 0.809 |
| Spiked # 6 | 0.660 |

Acknowledgements

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Chris Hudalla
- Jeremy Shia
- Kim van Tran