



# **Sugar & Carbohydrate Analysis**

## **A Practical Perspective**



# Liquid Chromatography System




# Sample Preparation




- Grind
- Extract
- Particles – Filtration
- Matrix elimination – extraction, dialysis, centrifuge, etc...




# Filtration



# Dialysis



# Dialysis



- particles
- high molecular weights
- high organic load
- oils
- fats
- proteins
- dyestuffs






# Carbohydrates: orange juice

Metrosep Carb1 – 250	ppm
1	n.d.
2	Inositol 1107
3	n.d.
4	n.d.
5	n.d.
6	Glucose 19214
7	Fructose 18820
8	Sucrose 26155

Dialysis, dilution 1:100

100 mmol/L NaOH




# Carbohydrates: milk

Metrosep Carb1 – 250	ppm
1-3 Glycerin	52
4 Inositol	35
5 Glycerin	28
6 Xylitol	21
7 Arabitol	6
8 Arabinose	51
9 Mannose	47
10 Glucose	53
11 Lactose	35625

Dilution 1:1000, dialysis

100 mmol/L NaOH



# Anions, pulsed amperometric detection

## Carbohydrates: milk chocolate




Metrosep Carb1 – 250

ppm		
3	Inositol	1.50
6	Arabitol	0.26
7	Sorbitol	0.40
8	Mannitol	0.77
10	Arabinose	1.29
11	Glucose	1.37
12	Xylose	2.09
13	Fructose	1.60
14	Lactose	838.96

Dilution 1:10, dialysis

100 mmol/L NaOH



- **Sugars and Carbohydrates fuel biological processes**
- **Reduce temperature in autosampler to slow / prevent change in sample**
- **Peltier cooling or external cooling**




- Specific for column and separation
- Isocratic for simple sugars
- Gradient for complex mixtures
- Follow manufacturer recommendations

**Column:** METROSEP Carb 1  
4X250mm

**Eluent:** 2.5 mmol/L sodium  
hydroxide + 0.65 mmol/L sodium  
acetate in ultrapure water

**Post column reagent:** 300 mmol/L  
sodium hydroxide in ultra pure water



- **Choice of vendors**
- **Selectivity/separation needed**
- **Isocratic and gradient separations**

# Column Selection

<u>Sugar Class</u>	<u>Separation Mode</u>
	Ligand Exchange +GFC
<b>Monosacharides Disacharides</b>	Ligand Exchange +Reversed Phase
	Normal Phase
<b>Sacharides &amp; Sugar Alcohols</b>	Ligand Exchange +GFC
	Ligand Exchange + Normal Phase
<b>Oligosacharides</b>	Ligand Exchange +GFC
<b>Amino Sugars</b>	Ligand Exchange +GFC
	Normal Phase
<b>Acidic Sugars</b>	Ion Exclusion
	Ion Exchange
<b>Sacharides &amp; Organic Acids</b>	Ion Exclusion+GFC
<b>Sugar Alcohols &amp; Organic Acids</b>	Ion Exclusion+GFC



# Isocratic Separation


**Sample:** Vodka

**Column:** Metrosep Carb 1  
4X150mm

**Mobile Phase:** 100mM NaOH

**Loop:** 20uL

**Detector:** Metrohm PAD





# Separation


**Sample:** Sugar Standards

**Column:** Metrosep Carb 1  
4x250mm

**Mobile Phase:** 200mM  
NaOH

**Loop:** 20uL

**Detector:** Metrohm PAD





**VARIAN**

## APPLICATION NOTE

### Organic acids and sugars


Application 1186 LC

#### Fermentation standard

Technique : HPLC  
Column : Metacarb 67H , 300 x 6.5 mm  
Part. No. A6244  
Mobile Phase : 0.025 N H<sub>2</sub>SO<sub>4</sub>  
Flow Rate : 0.6 mL/min  
Temperature : 65 °C  
Detection : RI  
Sample : 20 µL

#### Peak Identification

- 1 Maltotriose
- 2 Maltose
- 3 Glucose
- 4 Fructose
- 5 Lactic acid
- 6 Glycerol
- 7 Acetic acid
- 8 Methanol
- 9 Ethanol




# Separation



## Sugaralcohols, Mono- and Disaccharides with BioQuant Carbohydrate

Column: Sugar1  
Dimension: 250 x 4.6 mm PEEK  
Eluent: 100 mM NaOH  
Flow: 0.9 ml/min  
Detection: BioQuant Carbohydrate  
Temperature: 28 °C  
Injection: 20 µl  
Sample: 1 Inositol, 2 Arabitol, 3 Sorbitol, 4 Fucose, 5 Arabinose, 6 Glucose, 7 Xylose, 8 Fructose, 9 Lactose, 10 Sucrose (1-50 mg/l)




no. 171



## Carbohydrates, Organic Acids and Alcohols in Wine

### Carbohydrate H<sup>+</sup>

Part Number: 00253776  
Dimension: 300 x 7.8 mm  
Eluent: 1.25 mM H<sub>2</sub>SO<sub>4</sub>  
Flow: 0.6 ml/min  
Detection: RI  
Temperature: 45 °C  
Injection: 50 µl  
Sample: Wine Complet Standard





- Refractive Index
- Pulsed Amperometric
- Mass Spec



	<b>Selectivity</b>	<b>Sensitivity</b>	<b>Cost</b>	<b>Separation</b>
<b>Refractive Index</b>	Little	0.1 to 0.05%	Moderate	Isocratic
<b>Pulsed Amperometric</b>	Some	ng/g, ppb	Moderate	Isocratic, Gradient
<b>Mass Spec</b>	Most	ng/g, ppb	Expensive	Isocratic, Gradient

## Measuring principle

A defined decent potential is applied to a three electrode system consisting of reference-, working- and auxiliary electrodes. If any oxidizable substance pass these electrodes an oxidation process takes place. A current corresponding to the concentration of the respective substance is measured.

- Simple connection to any IC- or HPLC-system
- Large variety of electrode materials (glassy carbon, carbon paste, silver, gold, etc.)
- Highly sensitive
- Very selective


- **Electrospray ionization**
  - Operating temp 250 to 300°C
  - Sugars oxidize
  - Form adduct with HCl
- **HCl Adduct formation**
  - HCl is volatile
  - Source of  $\text{Cl}^-$  for adduction
  - Post column addition

## Chloride Adduction

- Carbohydrates chloride adduct and only these adducts are detectable in negative mode ESI.
- Chloride adducts of carbohydrates are collisionally stabilized in the API source.
- In carbohydrates, the chlorine appears to hydrogen bond to hydroxyl groups.
- Chloride adduction also occurs with alcohols such as glycerol, but the signal is weaker.
- Chloride does not adduct to organic acids and does not enhance the MS signal.



## Optimizing MS Response with Ramped Fragmentor Voltage



- Fragmentor voltage effects both sensitivity and fragmentation and is compound class specific
- Fragmentor voltage can change as a function of the mass (ramped) during scan or selected ion monitoring
- Ramped fragmentor voltage can be used to optimize the response from both molecular and fragment ions

- You have choices in hardware and separation
- Choose the right tools for your application
- What do you need?

# Thank You

**Larry Tucker**  
**Brinkmann Instruments, Inc**  
**ltucker@brinkmann.com**  
**281 330 9532**