

VICAM[®]

A Waters Business



New Approaches to Mycotoxin Control Programs: Simultaneous Rapid Detection of Several Mycotoxins

**2020 AAFCO Mid-year Meeting
Jan 22, 2020**

**Nancy Collette
Technical Service and Applications Manager
VICAM, A Waters Business**

- Mycotoxins frequently contaminate animal feeds.
- More than 70% of feeds analyzed globally are contaminated with mycotoxins. 38% with more than one mycotoxin.
- The monitoring of mycotoxins in feed ingredients can be used as a valuable tool for a quality control program for mycotoxins and protection of animal health.

OVERVIEW

1. Traditionally, mycotoxins are determined and regulated at the level of **feed ingredients**, facing difficulties such as the lack of homogeneity in the contamination of the lot and the transfer of the problem into feed. At this level, rapid methods for simultaneous multiple analyzes **in the field** are an invaluable tool.
2. Monitoring of **finished feeds** confirms toxin has not been transferred into the feed.
3. Biomonitoring of the **exposure of animals** to mycotoxins using the determination of risk biomarkers is a complementary tool in the control programs of mycotoxins.

Rapid and Simple monitoring of feed ingredients

With **one** short extraction step and **minimal** hands-on manipulation you can run **multiple** tests for **up to five** mycotoxins at once!



Rapid and Simple

Myco 5in1 - Single Extraction



Weigh 5g of finely-ground sample into an extraction tube, then add 25 mL of AQUA Premix



Extract the sample by vortexing at high speed for 2 minutes



Filter the extract into a clean tube

Afla- γ

Ochra- γ

Fumo- γ

Zearala- γ

DON- γ



Pipette 100 μ L of extract into the sample well at 1 drop/second



Place Don-V strip into mini-incubator at 40°C and close the lid for 2-10 minutes.



Pipette 100 μ L of extract into the sample well at 1 drop/second with the strip still inside the incubator



Incubate the strip for 5 minutes at 40°C



During the 5 minutes calibrate the Vertu by scanning the appropriate barcode




Place the strip inside the Vertu and press the center key to read




Place the strip on a flat surface for 5 minutes

Mycotoxin Level Low



Mycotoxin Level High



VICAM 5-in-1 Lateral Flow for Mycotoxin Diagnostics

Mycotoxin	LOD	Range
Aflatoxin	2 ppb	0-40 ppb, (<i>High range method to 100 ppb</i>)
Deoxynivalenol	0.25 ppm	0-4 ppm (1 st strip), to 16 ppm (2 nd strip)
Fumonsin	0.2 ppm	0-20 ppm (1 st strip), to 80 ppm (2 nd strip)
Ochratoxin	2 ppb	0-30 ppb
Zearalenone	0.1 ppm (100 ppb)	0-5 ppm

- Rapid: from sample to results in less than 10 minutes
- Easy to use- anyone can be trained to run the test
- Inexpensive
- Excellent for screening at the point of collection
- Meets worldwide regulatory limits.
- Non toxic - Water based extraction solution



USDA-FGIS Certifications

- Afla-V AQUA
- Fumo-V AQUA
 - DON-V

Real Results: Naturally Contaminated Corn

HPLC Results		Vertu Results					
		R1	R2	R3	Mean	SDr	RSDr%
Afla-V (ppb)	9.45	10	9.39	9.39	9.46	0.12	1.2
	18.3	17.1	18.9	18.9	18.30	1.0	5.6
DON-V (ppm)	1.8	1.62	1.82	1.80	1.75	0.11	6.28
	4.3	4.65	4.12	3.84	4.20	0.41	9.77
Fumo-V (ppm)	0.50	0.60	0.52	0.43	0.52	0.09	16.4
	4.8	4.08	4.56	4.52	4.39	0.27	6.0
Ochra-V (ppb)	4.2	5.48	4.46	5.81	5.25	0.70	13.4
	18.0	21.2	19.1	21.7	20.6	1.3	6.6
Zearala-V (ppm)	0.42	0.44	0.39	0.43	0.42	0.03	6.3
	1.36	1.42	1.51	1.45	1.46	0.05	3.1


- Mycotoxin testing **in the field** is best done by rapid quick in field assays such as **lateral flow strip tests on incoming ingredients such as corn**. Incoming grains can be rapidly screened before being made into a finished feed.
- Confirmation of mycotoxin results and testing of **finished feeds** and complex feed ingredients is best done by **laboratory instrumental methods such as HPLC, UPLC or LC/MS/MS**

Advantages of Immunoaffinity Chromatography clean up for HPLC, UPLC and LC/MS/MS for multiple mycotoxin analysis




- **Concentrates sample for lower limits of detection.** Other methods make sample more dilute (dilute and shoot) or do not concentrate sample.
- Gives more **toxin specific clean up for complex samples** such as feeds.
- **Reduces matrix materials** so reduces background peaks and interferences.
- Reduces amount of **cleaning needed for LC/MS/MS cone.**
- Also **reduces matrix enhancement or suppression on LC/MS/MS.**

Inside an immunoaffinity column



LC-MS (Liquid chromatography- Mass spectrometry)





Advantage of MS- It can detect molecules that do not fluoresce or absorb



Waters UPLC conditions – multi-toxin LC-MS/MS

- Waters ACQUITY UPLC I-Class with FL injector operating in PLNO mode
- Needle: PEEK, 10 μ L
- Column: ACQUITY UPLC BEH-C₁₈ column (2.1 \times 100mm, 1.7 μ m)
- Mobile phase A: methanol +0.5% acetic acid +0.1% formic acid
- Mobile phase B: 1mM ammonium acetate in water +0.5% acetic acid +0.1% formic
- Weak needle wash: H₂O:MeCN 1:1 +0.125mM EDTA Strong needle wash: H₂O +20mM citric acid:MeOH:MeCN : IPA: acetone: DMSO 37:9:19:19:9:7
- Flow rate = 0.4 mL min⁻¹
- Column temperature: 40° C




MS/MS conditions – multi-toxin method

Acquity UPLC®
Xevo™ TQ-XS

Parameter	Setting
Ionisation mode	ESI ^{+/-}
Capillary voltage kV	+0.5/-0.3
Source offset V	30
Source temperature °C	150
Desolvation temperature °C	500
Desolvation gas flow L/h	800
Cone gas flow L/h	150
Collision gas flow (mL/min) argon	0.15

2 transitions were monitored for each analyte & 1 transition was used for the respective ¹³C-labelled isomer



Characteristics of Mass Spectrometry

- Instrument can be used for other analysis besides mycotoxins (such as pesticides)
- Mass Spec can detect molecules that do not fluoresce or absorb
- Good for multiple mycotoxin analysis
- Good for confirmation of mycotoxins
- Must have a laboratory environment
- Triple Quad (MS/MS) gives best results (lower limits of detection and confirmation)
- Best for well trained scientists
- Need to accurately measure and store mycotoxins in order to get accurate results
- Need to use matrix matched calibration standards or internal calibrators to adjust for matrix affects

- 1. Lattanzio, V., Solfrizzo, M., Powers, S., Visconti, A., Simultaneous determination of aflatoxins, ochratoxin A and *Fusarium* toxins in maize by liquid chromatography/tandem mass spectrometry after multitoxin immunoaffinity cleanup, *Rapid Commun. Mass Spectrometry*, 2007; 21: 3253-3261.
- 2. Lattanzio, V., Ciasca, B., Powers, S., Visconti, A., Improved method for the simultaneous determination of aflatoxins, ochratoxin A and *Fusarium* toxins in cereals and derived products by liquid chromatography –tandem mass spectrometry after multi-toxin immunoaffinity clean up, *Journal of Chromatography A*, 2014; 1354: 139-143.
- 3. Park, J et al, Distribution Analysis of twelve Mycotoxins in Corn and Corn-Derived Products by LC-MS/MS to Evaluate the Carry-Over Ratio During Wet-Milling, *Toxins*, 2018; 10, 319.
- 4. Soleimany, F., Jinap, S., Rahmani, A., Khatib A., Simultaneous detection of 12 mycotoxins in cereals using RP-HPLC-PDA-FLD with PHRED and a post-column derivatization system, *Food Additives and Contaminants Part A Chem Anal Control Expo Risk Assess*, 2011 Apr 28(4):494-501.
- 5. Tang, Y.Y., Lin H.Y., Chen Y.C., Su, W.T. Wang. S.C., Chiueh L.C., Shin Y.C., Development of a Quantitative Multi-Mycotoxin Method in Rice, Maize, Wheat and Peanut Using UPLC-MS/MS, *Food Anal. Methods*, 2013; 6:727-736.
- 6. Kim, D.H., Jang, H.S., Choi, G.H., Kim, H.J., Kim, H.J., Kim, H.L., Cho, H.J., Lee C., Occurrence of Mycotoxins in Korean Grains and Their Simultaneous Analysis, *Korean J. Food Sci. Technol.*, 2013; Vol. 45, No 1, pp. 111-119.
- 7. Solfrizzo, M., Gambacorta, L., Lattanzio V.M.T., Powers, S., Visconti, A., Simultaneous LC–MS/MS determination of aflatoxin M1, ochratoxin A, deoxynivalenol, de-epoxydeoxynivalenol, α and β -zearalenols and fumonisin B1 in urine as a multi-biomarker method to assess exposure to mycotoxins, *Anal. Bioanal. Chem.*, 2011; 401:2831-2841.
- 8. Vaclavikova, M., MacMahon, S., Zhang, K., Begley, T. Application of single immunoaffinity clean-up for simultaneous determination of regulated mycotoxins in cereals and nuts, *Talanta*, 2013; 117: 345-351.
- 9. Kim, D.H., et al, Simultaneous Determination of Multi-Mycotoxins in Cereal Grains Collected from South Korea by LC/MS/MS, *Toxins*, 2017; 9, 106 : 1-13.

In addition to testing for mycotoxins in incoming ingredients and finished feed, **exposure** of dairy cattle to aflatoxin can be determined by testing for aflatoxin M₁ in their milk.

Monitoring exposure of dairy cows to aflatoxins

- Aflatoxin M₁ is the metabolic byproduct created when a dairy animal ingests feed that contains aflatoxin B₁.
- Aflatoxin M₁ is a Group 2B carcinogen.
- Aflatoxin M₁ is regulated worldwide in raw and processed milk:



Contaminated Feed Rations...



**Hydroxylation
of AFB₁/AFB₂**



**AFM₁:
Classified as
a Group 2B
Carcinogen
by the IARC**



**Got [AFM₁
Free] Milk ?**


**AFM₁
Predominant
Residue: Heat
Stable, Survives
Processing**



**Secreted as
Less Toxic
AFM₁/AFM₂
in Milk**



World AFM₁ Thresholds



BIOMARKER

Metabolite resulting from transformation of a molecule to which an individual has been exposed and which can be identified in body fluids or excreta of the receiving organism.

Dairy Herd Improvement



AFLA B₁/B₂ in Maize,
Cottonseed, CGM,
Distillers Grains,
Hominy Feed, Peanut
Meal, Canola, Silage,
Copra, Citrus Pulp,
Soybean Meal, Wheat
Bran, Rice



Decreased Milk
Production, Ataxia,
Elevated SCCs',
Feed Refusal &
Reduced
Reproductive
Efficiency



Afla M₁-V
QUANTITATIVE STRIP TESTS

Opportunity for Aflatoxin M₁



← AFB₁ → AFM₁: 1 to 6% Carry-Over to Milk, 12 - 24 hours from first ingestion →

Aflatoxin M₁ as a Biomarker for aflatoxin B₁ exposure

- Aflatoxin M₁ appears in milk at the first milking (1hr) after cows are fed aflatoxin B₁ contaminated feed
- After feeding aflatoxin B₁ ceases, aflatoxin M₁ decreases in the milk for up to 7 days
- Testing for aflatoxin M₁ in milk can determine if cows have been exposed to aflatoxin B₁ in their feed



Monitoring exposure of cows to aflatoxins

- Biomarkers can not only be indicators of mycotoxin exposure, but also excellent and invaluable indicators of the effectiveness of binders added to feed.
- The methodologies for the determination of biomarkers are especially valuable if they can be used in the field in a non-invasive, fast and simple way.

Rapid and Simple




**Add 200µl
Cold Milk**

**Vortex 3
times for 5
seconds
each time
to mix**

**Add strip
and
incubate
40°C for 10
minutes**

**Read
on Vertu**

**Limit of Detection 25ppt (0.025ppb)
Assay range 0.025-0.75ppb**



Rapid Milk AFM₁ Testing

Lateral Flow



- Dairy Farm
- Transfer Stations
- Receiving Stations
- Milk Depots
- Milk Collection Centers

Lateral Flow



- Milk Receiving Bay – Processing Plant
- Milk Testing Lab – Processing Plant





Lateral Flow, ELISA, Fluorometric-IAC, LC-IAC





- Central Milk Testing/ Producer Payment Lab
- Import Inspection Lab
- Regulatory / Public Health Lab
- Dairy Extension/ Research Lab
- Corporate Industry Lab, Contract Testing Lab
- Bovine Veterinary Producer Support Lab

Test Performance – Linearity & Range

Whole Homogenized Cow Milk




Raw Cow Milk




Range covers
most worldwide
regulations for
Afla M₁ levels in
one rapid test

Test Performance – Linearity & Range


Pasteurized Skimmed Cow Milk




Pasteurized Whole Cow Milk (non-homogenized)




Raw Goat Milk



Raw Sheep Milk



Test Performance – Precision and Repeatability



Afla M₁-V kit part # 176004148



Afla M₁-V BEQ Part # 176004172 (110v) 176004173 (220V)



Conclusions


- The combined use of methodologies for multiple simultaneous analyzes of mycotoxins in feed ingredients and analysis of biomarkers in fluids of animals that consume the feed, may in the immediate future improve the control of problems arising from exposure to mycotoxins.
- Would you like information on multiply mycotoxin testing in corn or feed or aflatoxin M₁ testing in milk?
- Please contact me here in person or after the meeting at nancy_zabe_collette@waters.com

Any
Questions?



Thank You

When mycotoxin contamination is positive



When sample is mycotoxin free (negative)

