

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
300.01	Aflatoxin, Neogen Vera-Tox	PPB		
300.02	Aflatoxin, Neogen Agri-Screen	PPB	990.32	49.2.16
300.03	Aflatoxin, Vicam Aflatest.....	PPB	991.31	49.2.18
300.04	Aflatoxin, Rhone Poulenc Aflascan	PPB		
300.05	Aflatoxin, Charm	PPB		
300.06	Aflatoxin, Romer Minicolumn	PPB	975.36	49.2.05
300.07	Aflatoxin, Biopharm Ridascreen	PPB		
300.10	Aflatoxin, TLC	PPB		
300.20	Aflatoxin, HPLC	PPB		
300.99	Aflatoxin, Misc.	PPB		
120.00	Alanine, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
120.01	Alanine, Pre-col OPA Der.....	PCT		
120.02	Alanine, Post-col OPA Der	PCT		
120.03	Alanine, Pre-col Dansyl Cl Der	PCT		
120.04	Alanine, Pre-col FMOC Der	PCT		
120.05	Alanine, Pre-col AQC Der.....	PCT		
120.99	Alanine, Misc.	PCT		
015.00	Aluminum, ICP	PPM		
015.99	Aluminum, Misc.	PPM		
045.00	Amprolium, Color.....	PCT	961.24	5.1.08
045.01	Amprolium, Fluor	PCT	965.47	5.1.09
045.02	Amprolium, HPLC.....	PCT		
045.99	Amprolium, Misc.	PCT		
046.00	Apramycin, Micro	g/Ton		
046.02	Apramycin, HPLC.....	g/Ton		
046.99	Apramycin, Misc.	g/Ton		
200.01	Arachidonic Acid.....	PCT		
121.00	Arginine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
121.01	Arginine, Pre-col OPA Der.....	PCT		
121.02	Arginine, Post-col OPA Der	PCT		
121.03	Arginine, Pre-col Dansyl Cl Der	PCT		
121.04	Arginine, Pre-col FMOC Der.....	PCT		
121.05	Arginine, Pre-col AQC Der	PCT		
121.99	Arginine, Misc.	PCT		
047.00	Arsanilic Acid.....	PCT	954.17	5.1.12
047.01	Arsanilic Acid, Total As	PCT	957.22	5.1.13
047.02	Arsanilic Acid, HPLC	PCT		
047.99	Arsanilic Acid, Misc.....	PCT		
016.00	Arsenic, AA, Hydride	PPM		
016.01	Arsenic, Total	PPM	957.22	5.1.13
016.02	Arsenic, ICP	PPM		
016.99	Arsenic, Misc.	PPM		
005.00	Ash	PCT	942.05	4.1.10
005.01	Ash, Sugars & Syrups	PCT	900.02	44.1.05
005.02	Ash, LECO	PCT		
005.03	Ash, Microwave furnace	PCT		
005.11	Ash, NIR	PCT		
005.99	Ash, Misc.	PCT		
122.00	Aspartic, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
122.01	Aspartic, Pre-col OPA Der.....	PCT		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
122.02	Aspartic, Post-col OPA Der	PCT		
122.03	Aspartic, Pre-col Dansyl Cl Der	PCT		
122.04	Aspartic, Pre-col FMOC Der	PCT		
122.05	Aspartic, Pre-col AQC Der	PCT		
122.99	Aspartic, Misc.	PCT		
048.00	Bacitracin, Plate	g/Ton	957.24	5.3.03
048.01	Bacitracin, Plate MeOH Ext	g/Ton	993.29	5.3.05
048.02	Bacitracin, Turbid	g/Ton		
048.04	Bacitracin, HPLC	g/Ton	982.44	5.2.01
048.99	Bacitracin, Misc.	g/Ton		
049.00	Bambermycin, Plate	g/Ton		
049.99	Bambermycin, Misc.	g/Ton		
040.00	Barium, ICP	PPM		
040.99	Barium, Misc.	PPM		
114.01	Biotin, Microbiological	mg/KG		
017.00	Boron, ICP	PPM		
017.99	Boron, Misc.	PPM		
018.01	Cadmium, AA	PPM		
018.02	Cadmium, ICP	PPM		
018.99	Cadmium, Misc.	PPM		
019.00	Calcium, Ox-Mn04 Vol.	PCT	927.02	4.8.03
019.01	Calcium, AA Spect	PCT	968.08	4.8.02
019.02	Calcium, Hach Method	PCT		
019.03	Calcium, Semiauto (Autoanalyzer)	PCT		
019.04	Calcium, Em Spect	PCT		
019.05	Calcium, ICP, Dry Ash	PCT		
019.06	Calcium, DCP	PCT		
019.07	Calcium, PreWet	PCT		
019.08	Calcium, EDTA	PCT		
019.09	Calcium, ICP, Wet Ash	PCT		
019.99	Calcium, Misc.	PCT		
050.00	Carbadox, Color	PCT	977.35	5.1.17
050.01	Carbadox, HPLC	PCT		
050.99	Carbadox, Misc.	PCT		
100.00	Carotene	mg/LB	970.64	45.1.04
100.01	Carotene, HPLC	mg/LB		
100.99	Carotene, Misc.	mg/LB		
051.00	Chlortetracycline, Plate	g/Ton	967.39	5.3.06
051.01	Chlortetracycline, Turbid	g/Ton	977.37	5.3.07
051.02	Chlortetracycline, Chem Premix	g/Ton	JAOAC <u>47</u> & <u>46</u>	
051.03	Chlortetracycline, HPLC	g/Ton		
051.99	Chlortetracycline, Misc.	g/Ton		
101.00	Choline Chloride, Micro	mg/LB		
101.01	Choline Chloride, Chem	mg/LB		
101.02	Choline Chloride, HPLC	mg/LB		
101.99	Choline Chloride, Misc.	mg/LB		
020.00	Chromium, AA	PPM		
020.01	Chromium, ICP	PPM		
020.99	Chromium, Misc.	PPM		
052.00	Clopidol	g/Ton		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
052.99	Clopidol, Misc.	g/Ton		
021.00	Cobalt, Color	PPM	952.02	4.8.07
021.01	Cobalt, AA	PPM	968.08	4.8.02
021.02	Cobalt, ICP	PPM		
021.03	Cobalt, DCP	PPM		
021.99	Cobalt, Misc.	PPM		
022.00	Copper, Color	PPM	947.03	4.8.08
022.01	Copper, AA	PPM	968.08	4.8.02
022.02	Copper, Em Spect	PPM	953.01	3.2.01
022.03	Copper, ICP, Dry Ash	PPM	968.08 (D)(a)	
022.04	Copper, DCP	PPM		
022.05	Copper, ICP, Wet Ash	PPM	968.08 (D)(b)	
022.99	Copper, Misc.	PPM		
053.00	Coumaphos, GLC	PCT		
053.01	Coumaphos, Polarographic	PCT		
053.99	Coumaphos, Misc.	PCT		
124.00	Cysteine/Cystine, PAO Post-col Ninhydri	PCT	994.12	4.1.11
124.01	Cysteine/Cystine, PAO Pre-col OPA Der	PCT		
124.02	Cysteine/Cystine, PAO Post-col OPA Der	PCT		
124.03	Cysteine/Cystine, PAO Pre-col Dansyl Cl	PCT		
124.04	Cysteine/Cystine, PAO Pre-col FMOC Der	PCT		
124.05	Cysteine/Cystine, PAO Pre-col AQC Der	PCT		
124.99	Cysteine/Cystine, Misc.	PCT		
054.00	Decoquinatate	mg/LB	969.55	5.1.18
054.01	Decoquinatate, HPLC	mg/LB		
054.99	Decoquinatate, Misc.	mg/LB		
055.00	Erythromycin, Plate	g/Ton	971.48	5.3.08
055.99	Erythromycin, Misc.	g/Ton		
056.00	Ethopabate	PCT	964.29	5.1.22
056.01	Ethopabate, HPLC UV	PCT		
056.02	Ethopabate, HPLC, Fluor	PCT		
056.99	Ethopabate, Misc.	PCT		
057.00	Ethoxyquin	PCT	963.07	4.10.02
057.99	Ethoxyquin, Misc.	PCT		
003.00	Fat, Eth Ext., Direct	PCT	920.39	4.5.01
003.01	Fat, Ind Eth Ext (13 th ed.), Indirect	PCT	920.39	4.5.01
003.04	Fat, In Fish Meal	PCT	948.04	4.5.04
003.05	Fat, Distillers Meth	PCT		
003.06	Fat, Pet Ether	PCT		
003.07	Fat, Aq Ext	PCT	920.39	4.5.01
013.08	Fat, Roese-Gottlieb Modified	PCT		
003.09	Fat, Soxtec, Eth Ext	PCT	2003.05	4.5.05
003.10	Fat, Soxtec, Pet Ether	PCT		
003.11	Fat, NIR	PCT		
003.12	Fat, Hexane Ext	PCT		
003.13	Fat, Soxtec, Hexane Ext	PCT	2003.06	4.5.06
003.14	Fat, Ankom	PCT		
003.99	Fat, Misc.	PCT		
013.02	Fat, Mojonier, Bak Ext	PCT	954.02	4.5.02
013.03	Fat, Roese-Gottlieb	PCT	932.02	4.5.03

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
013.10	Fat, Soxtec-Acid Hydrolysis.....	PCT		
013.11	Fat, Super Critical Fluid Extraction	PCT		
013.12	Fat, NIR-Acid Hydrolysis	PCT		
013.13	Fat, Ankom-Acid Hydrolysis	PCT		
013.99	Fat, pretreat or extended ext., Misc.	PCT		
092.00	Fenbendazole.....	PPM		
004.00	Fiber, Crude Asbestos Free	PCT	962.09	4.6.01
004.01	Fiber, Sing Filt	PCT		
004.03	Fiber, Fritted Glass	PCT	978.10	4.6.02
004.06	Fiber, Fibertec	PCT		
004.07	Fiber, ANKOM	PCT		
004.11	Fiber, NIR	PCT		
004.99	Fiber, Misc.	PCT		
008.02	Fiber, Acid Detergent.....	PCT	973.18	4.6.03
008.05	Fiber, Acid Detergent-Hach	PCT		
008.08	Fiber, Acid Detergent ANKOM	PCT		
008.99	Fiber, Acid Detergent Misc.	PCT		
009.04	Fiber, Neutral Det-No ENZ Pretreat	PCT		
009.07	Fiber, Neutral Det-Amylase Pretreat	PCT	2002.04	4.6.04
009.09	Fiber, Neutral Det ANKOM	PCT		
009.99	Fiber, Neutral Det Misc.	PCT		
023.00	Fluorine, Color	PCT	944.08	9.2.11
023.01	Fluorine, Ion Sel Electrode	PCT	975.08	4.8.09
023.99	Fluorine, Misc.	PCT		
113.01	Folic Acid.....	mg/KG	944.12	45.2.03
058.00	Furazolidone.....	PCT	960.63	5.1.24
058.01	Furazolidone, HPLC	PCT	985.51	5.1.23
058.99	Furazolidone, Misc.	PCT		
123.00	Glutamic Acid	PCT		
123.99	Glutamic Acid, Misc.	PCT		
125.00	Glutamic, Post-col Ninhydrin Der	PCT	994.12	4.1.11
125.01	Glutamic, Pre-col OPA Der.....	PCT		
125.02	Glutamic, Post-col OPA Der	PCT		
125.03	Glutamic, Pre-col Dansyl Cl Der.....	PCT		
125.04	Glutamic, Pre-col FMOC Der.....	PCT		
125.05	Glutamic, Pre-col AQC Der.....	PCT		
125.99	Glutamic, Misc.	PCT		
126.00	Glycine, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
126.01	Glycine, Pre-col OPA Der.....	PCT		
126.02	Glycine, Post-col OPA Der	PCT		
126.03	Glycine, Pre-col Dansyl Cl Der	PCT		
126.04	Glycine, Pre-col FMOC Der.....	PCT		
126.05	Glycine, Pre-col AQC Der.....	PCT		
126.99	Glycine, Misc.	PCT		
127.00	Histidine, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
127.01	Histidine, Pre-col OPA Der	PCT		
127.02	Histidine, Post-col OPA Der	PCT		
127.03	Histidine, Pre-col Dansyl Cl Der	PCT		
127.04	Histidine, Pre-col FMOC Der	PCT		
127.05	Histidine, Pre-col AQC Der.....	PCT		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
127.99	Histidine, Misc.	PCT		
059.00	Hygromycin B	UNITS/LB	960.67	5.3.09
059.99	Hygromycin B, Misc.	UNITS/LB		
024.00	Iodine, Knaph-Lamb	PPM	934.02	4.8.10
024.01	Iodine, Elm-Cald	PPM	935.14	4.8.11
024.02	Iodine, Iodized Salt	PPM	925.56	11.2.02
024.03	Iodine, Ion Sel Electrode	PPM		
024.99	Iodine, Misc.	PPM		
060.00	Ipronidazole	PCT	971.44	5.1.27
060.01	Ipronidazole, HPLC	PCT		
060.99	Ipronidazole, Misc.	PCT		
025.00	Iron, Color	PPM	935.12	4.8.01
025.01	Iron, AA	PPM	968.08	4.8.02
025.02	Iron, Em Spect	PPM	953.01	3.2.01
025.03	Iron, ICP, Dry Ash	PPM	968.08 (D)(a)	
025.04	Iron, DCP	PPM		
025.05	Iron, ICP, Wet Ash	PPM	968.08 (D)(b)	
025.99	Iron, Misc.	PPM		
128.00	Isoleucine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
128.01	Isoleucine, Pre-col OPA Der	PCT		
128.02	Isoleucine, Post-col OPA Der	PCT		
128.03	Isoleucine, Pre-col Dansyl Cl Der	PCT		
128.04	Isoleucine, Pre-col FMOC Der	PCT		
128.05	Isoleucine, Pre-col AQC Der	PCT		
128.99	Isoleucine, Misc.	PCT		
130.00	L-Lysine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
130.01	L-Lysine, Pre-col OPA Der	PCT		
130.02	L-Lysine, Post-col OPA Der	PCT		
130.03	L-Lysine, Pre-col Dansyl Cl Der	PCT		
130.04	L-Lysine, Pre-col FMOC Der	PCT		
130.05	L-Lysine, Pre-col AQC Der	PCT		
130.99	L-Lysine, Misc.	PCT		
061.00	Lasalocid, Sodium (Microbio)	g/Ton	975.60	5.3.10
061.01	Lasalocid, Sodium (Spect)	g/Ton	975.61	5.2.04
061.02	Lasalocid, Sodium (HPLC)	g/Ton		
061.99	Lasalocid, Sodium, Misc.	g/Ton		
026.00	Lead	PPM		
026.99	Lead, Misc.	PPM		
129.00	Leucine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
129.01	Leucine, Pre-col OPA Der	PCT		
129.02	Leucine, Post-col OPA Der	PCT		
129.03	Leucine, Pre-col Dansyl Cl Der	PCT		
129.04	Leucine, Pre-col FMOC Der	PCT		
129.05	Leucine, Pre-col AQC Der	PCT		
129.99	Leucine, Misc.	PCT		
062.00	Levamisole, HCl, GLC	PCT		
062.01	Levamisole, HCl, HPLC	PCT		
062.99	Levamisole, Misc.	PCT		
063.00	Lincomycin, Plate	g/Ton	967.40	5.3.11
063.99	Lincomycin, Misc.	g/Ton		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
210.01	Linoleic Acid	PCT		
001.00	Loss on Drying, Vac 95 deg 5 hr	PCT	934.01	4.1.03
001.02	Loss on Drying, Vac on sand	PCT	925.45	44.1.03D
001.03	Loss on Drying, ISO 6496	PCT		
001.05	Loss on Drying, LECO	PCT		
001.06	Loss on Drying, Vac on pumice stone	PCT	925.45	44.1.03C
001.07	Loss on Drying, 104 deg 3 hr., In malt	PCT	935.29	27.3.06
001.08	Loss on Drying, 102 deg 16 hr., In meat	PCT	950.46	39.1.02B
001.99	Loss on Drying, Misc.	PCT		
011.01	Loss on Drying, 135 deg 2 hr	PCT	930.15	4.1.06
011.99	Loss on Drying, High Temp. Methods, Misc.	PCT		
140.00	Lysine, Free (Available).....	PCT	975.44	45.4.03
027.00	Magnesium, Grav MG2P207	PCT	937.01	2.6.19
027.01	Magnesium, AA	PCT	968.08	4.8.02
027.02	Magnesium, Em Spect	PCT	953.01	3.2.01
027.03	Magnesium, ICP, Dry Ash	PCT	968.08 (D)(a)	
027.04	Magnesium, DCP	PCT		
027.05	Magnesium, ICP, Wet Ash	PCT	968.08 (D)(b)	
027.06	Magnesium, Hach Method	PCT		
027.99	Magnesium, Misc.	PCT		
028.00	Manganese, Color	PPM	917.04	4.8.12
028.01	Manganese, AA	PPM	968.08	4.8.02
028.02	Manganese, Em Spect	PPM		
028.03	Manganese, ICP, Dry Ash	PPM	968.08 (D)(a)	
028.04	Manganese, DCP	PPM		
028.05	Manganese, ICP, Wet Ash	PPM	968.08 (D)(b)	
028.99	Manganese, Misc.	PPM		
064.00	Melengestrol Acetate, GLC	PCT	980.36	5.1.28
064.01	Melengestrol Acetate, HPLC	PCT		
064.99	Melengestrol Acetate, Misc.	PCT		
029.00	Mercury.....	PPM		
029.99	Mercury, Misc.	PPM		
131.00	Methionine, PAO Post-col Ninhydrin Der	PCT	994.12	4.1.11
131.01	Methionine, PAO Pre-col OPA Der.....	PCT		
131.02	Methionine, PAO Post-col OPA Der	PCT		
131.03	Methionine, PAO Pre-col Dansyl Cl.....	PCT		
131.04	Methionine, PAO Pre-col FMOC Der.....	PCT		
131.05	Methionine, PAO Pre-col AQC Der	PCT		
131.99	Methionine, Misc.	PCT		
010.03	Moisture, Karl-Fischer	PCT	2001.12	4.1.09A
010.04	Moisture, Toluene Dist.....	PCT	925.04	4.1.04
010.11	Moisture, NIR.....	PCT		
010.99	Moisture, Misc.	PCT		
038.00	Molybdenum, ICP	PPM		
038.99	Molybdenum, Misc.	PPM		
065.00	Monensin, Plate.....	g/Ton	972.56	5.3.13
065.01	Monensin, Turbid.....	g/Ton	976.37	5.3.14
065.03	Monensin, HPLC	g/Ton	997.04	5.2.05
065.99	Monensin, Misc.	g/Ton		
000.01	NPN, Urea + Am, Urease Method	PCT	941.04	4.3.01

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
000.03	NPN, Automated.....	PCT		
091.00	Narasin	PCT		
066.00	Neomycin.....	g/Ton	970.89	5.3.15
066.99	Neomycin, Misc.	g/Ton		
102.00	Niacin, Chem	mg/LB	961.14	45.1.10
102.01	Niacin, Micro.....	mg/LB	944.13	45.2.04
102.02	Niacin, HPLC	PCT		
102.99	Niacin, Misc.	PCT		
067.00	Nicarbazin.....	PCT	956.11	5.1.30
067.99	Nicarbazin, Misc.	PCT		
039.01	Nickel, AA.....	PPM		
039.02	Nickel, ICP.....	PPM		
039.99	Nickel, Misc.	PPM		
030.00	Nitrate, Color	PCT	968.07	4.4.03
030.01	Nitrate, Ion Sel Electrode.....	PCT		
030.99	Nitrate, Misc.	PCT		
068.00	Nitrofurazone	PCT	960.63	5.1.24
068.01	Nitrofurazone, HPLC	PCT		
068.99	Nitrofurazone, Misc.	PCT		
069.00	Nitromide, 3,5-Dinitrobenzamide	PCT	959.18	5.1.37
069.01	Nitromide, HPLC	PCT		
069.99	Nitromide, Misc.	PCT		
070.00	Novobiocin.....	PCT	962.25	5.3.16
070.99	Novobiocin, Misc.	PCT		
071.00	Nystatin	PCT	974.48	5.3.17
071.99	Nystatin, Misc.	PCT		
072.00	Ormetroprim	PCT		
072.99	Ormetroprim, Misc.	PCT		
073.00	Oxytetracycline, Plate.....	g/Ton	968.50	5.3.19
073.01	Oxytetracycline, Turbid.....	g/Ton		
073.02	Oxytetracycline, High Mineral Feed.....	g/Ton		
073.03	Oxytetracycline, HPLC	g/Ton		
073.99	Oxytetracycline, Misc.	g/Ton		
103.00	Pantothenic Acid, Titrimetric	mg/LB	945.74	45.2.05
103.01	Pantothenic Acid, Microbiological.....	mg/LB	945.74	45.2.05
103.02	Pantothenic Acid, HPLC	mg/LB		
103.99	Pantothenic Acid, Misc.	mg/LB		
074.00	Penicillin, Plate	g/Ton	967.41	5.3.20
074.01	Penicillin, Ch Premix.....	g/Ton		
074.02	Penicillin, HPLC.....	g/Ton		
074.99	Penicillin, Misc.	g/Ton		
075.00	Phenothiazine	PCT	967.37	5.1.39
075.99	Phenothiazine, Misc.	PCT		
132.00	Phenylalanine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
132.01	Phenylalanine, Pre-col OPA Der	PCT		
132.02	Phenylalanine, Post-col OPA Der.....	PCT		
132.03	Phenylalanine, Pre-col Dansyl Cl Der	PCT		
132.04	Phenylalanine, Pre-col FMOC Der	PCT		
132.05	Phenylalanine, Pre-col AQC Der	PCT		
132.99	Phenylalanine, Misc.	PCT		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
031.00	Phosphorus, Vol	PCT	964.06	4.8.13
031.01	Phosphorus, Photometric	PCT	965.17	4.8.14
031.02	Phosphorus, GQMP (2.028)	PCT	964.06	4.8.13
031.03	Phosphorus, Autoanalyzer	PCT		
031.04	Phosphorus, Em Spect.....	PCT	953.01	3.2.01
031.05	Phosphorus, ICP	PCT		
031.06	Phosphorus, Hach Method	PCT		
031.07	Phosphorus, DCP	PCT		
031.08	Phosphorus, PreWet	PCT		
031.99	Phosphorus, Misc.	PCT		
150.00	Phytase, Colorimetric	FTU/LB	2000.12	4.10.06
150.99	Phytase, Misc.	FTU/LB		
076.00	Piperazine	PCT	967.38	5.1.40
076.99	Piperazine, Misc.	PCT		
032.00	Potassium	PCT	949.01	2.5.02
032.01	Potassium, AA	PCT	975.03	3.2.05
032.02	Potassium, Flame Emission	PCT	956.01	3.3.13
032.03	Potassium, STPB	PCT		
032.04	Potassium, Em Spect	PCT	953.01	3.2.01
032.05	Potassium, ICP	PCT		
032.06	Potassium, DCP	PCT		
032.07	Potassium, Hach Method	PCT		
032.99	Potassium, Misc.	PCT		
133.00	Proline, Post-col Ninhydrin Der	PCT	994.12	4.1.11
133.03	Proline, Pre-col Dansyl Cl Der	PCT		
133.04	Proline, Pre-col FMOC Der.....	PCT		
133.05	Proline, Pre-col AQC Der	PCT		
133.99	Proline, Misc.	PCT		
002.00	Protein, Crude	PCT	954.01	4.2.02
002.01	Protein, Auto Kjell-Foss	PCT	976.05	4.2.05
002.02	Protein, Semiauto Autoanalyzer	PCT	976.06	4.2.06
002.03	Protein, Hach Method	PCT		
002.04	Protein, Copper Cat	PCT	984.13	4.2.09
002.05	Protein, Copper, Boric Acid	PCT	2001.11	4.2.11
002.06	Protein, Combustion Nitrogen Analyzer	PCT	990.03	4.2.08
002.07	Protein, Block Dig	PCT	976.06	4.2.06
002.08	Protein, Cu/Ti	PCT	988.05	4.2.03
002.09	Protein, Selenium Catalyst	PCT		
002.10	Protein, Block dig/distillation	PCT		
002.11	Protein, NIR	PCT		
002.50	Protein, Crude Dig % Crude Protein.....	PCT	971.09	4.4.04
002.51	Protein, Pepsin Ind % Crude Protein.....	PCT	971.09	4.4.04
002.52	Protein, Pepsin Ind Residue	PCT	971.09	4.4.04
002.99	Protein, Misc.	PCT		
077.00	Pyrantel Tartrate	g/Ton	978.30	5.1.41
077.01	Pyrantel Tartrate, HPLC	g/Ton		
077.99	Pyrantel Tartrate, Misc.	g/Ton		
112.00	Pyridoxine (Vitamin B6)	mcg/G	961.15	45.2.08
112.99	Pyridoxine, Misc.	mcg/G		
093.00	Ractopamine Hydrochloride	g/Ton		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
104.00	Riboflavin, Fluorometric	mg/LB	970.65	45.1.08
104.01	Riboflavin, Micro	mg/LB	940.33	45.2.06
104.02	Riboflavin, Micro, Turbid.....	mg/LB	940.33	45.2.06
104.03	Riboflavin, HPLC	mg/LB		
104.99	Riboflavin, Misc.	mg/LB		
078.00	Roxarsone, Total As	g/Ton	957.22	5.1.13
078.01	Roxarsone, Color.....	g/Ton	971.47	5.1.46
078.02	Roxarsone, HPLC.....	g/Ton		
078.99	Roxarsone, Misc.	g/Ton		
079.00	Salinomycin, Plate.....	g/Ton		
079.01	Salinomycin, HPLC.....	g/Ton		
079.02	Salinomycin, Turbid	g/Ton		
079.99	Salinomycin, Misc.	g/Ton		
033.00	Salt, Sol Cl.....	PCT	943.01	4.8.05
033.01	Salt, Poten Cl	PCT	969.10	4.8.06
033.03	Salt, Quantab.....	PCT		
033.04	Salt, Hach Method	PCT		
033.05	Salt, Ion Sel Electrode	PCT		
033.99	Salt, Misc.	PCT		
034.00	Selenium, DCP	PPM		
034.01	Selenium, Fluor	PPM	969.06	3.4.13
034.02	Selenium, Grav, 11th.....	PPM	935.05	3.4.05
034.03	Selenium, AA, Flame.....	PPM		
034.04	Selenium, AA, Hydride	PPM		
034.05	Selenium, ICP	PPM		
034.06	Selenium, AA, Furnace.....	PPM		
034.99	Selenium, Misc.	PPM		
134.00	Serine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
134.01	Serine, Pre-col OPA Der	PCT		
134.02	Serine, Post-col OPA Der.....	PCT		
134.03	Serine, Pre-col Dansyl Cl Der.....	PCT		
134.04	Serine, Pre-col FMOC Der	PCT		
134.05	Serine, Pre-col AQC Der	PCT		
134.99	Serine, Misc.	PCT		
035.00	Sodium, AA.....	PCT		
035.01	Sodium, Ion Sel Electrode	PCT		
035.02	Sodium, Em Spect	PCT	953.01	3.2.01
035.03	Sodium, ICP	PCT		
035.04	Sodium, DCP.....	PCT		
035.05	Sodium, Flame Emission.....	PCT	956.01	3.3.13
035.99	Sodium, Misc.	PCT		
012.00	Starch, Polarimetric (Ewers).....	PCT		
012.01	Starch, Megazyme.....	PCT	996.11	32.2.05A
012.02	Starch, Colorimetric (GOP).....	PCT		
012.03	Starch, Enzymatic.....	PCT		
012.04	Starch, YSI Analyzer	PCT		
012.11	Starch, NIR.....	PCT		
012.99	Starch, Misc.....	PCT		
080.00	Streptomycin.....	g/Ton	971.49	5.3.22
080.99	Streptomycin, Misc.	g/Ton		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
006.00	Sugar, Sucrose	PCT	925.05	4.7.02
006.01	Sugar, Mod. Fehling Soln.	PCT	974.06	4.7.01
006.02	Sugar, TSI L-E C-V.....	PCT	925.05	4.7.02
006.03	Sugar, Invert w/o Inversn	PCT	970.58	44.2.09
006.04	Sugar, UNF Red Sub	PCT	948.23	44.2.10
006.05	Sugar, TSI, Lane-Eynon (12th).....	PCT	923.09	44.1.15
006.99	Sugar, Misc.	PCT		
081.00	Sulfadimethoxine	PCT		
081.99	Sulfadimethoxine, Misc.	PCT		
082.00	Sulfamethazine.....	PCT	969.57	5.1.50
082.01	Sulfamethazine, HPLC	PCT		
082.02	Sulfamethazine, HPLC-PCD	PCT	999.16	5.1.58
082.99	Sulfamethazine, Misc.	PCT		
083.00	Sulfanitran, APNPS	PCT	963.34	5.1.51
083.01	Sulfanitran, HPLC.....	PCT		
083.99	Sulfanitran, Misc.	PCT		
084.00	Sulfaquinoxaline	PCT	963.35	5.1.52
084.01	Sulfaquinoxaline, HPLC.....	PCT		
084.99	Sulfaquinoxaline, Misc.	PCT		
085.00	Sulfathiazole, Color	PCT		
085.01	Sulfathiazole, HPLC	PCT		
085.02	Sulfathiazole, HPLC-PCD.....	PCT	999.16	5.1.58
085.99	Sulfathiazole, Misc.	PCT		
036.00	Sulfur (Gravimetric)	PCT		
036.01	Sulfur (Analyzer).....	PCT		
036.02	Sulfur (Turbidimetric)	PCT		
036.03	Sulfur, ICP	PCT		
036.04	Sulfur, LECO	PCT		
036.99	Sulfur, Misc.	PCT		
139.00	Taurine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
139.01	Taurine, Pre-col OPA Der.....	PCT		
139.02	Taurine, Post-col OPA Der	PCT		
139.03	Taurine, Pre-col Dansyl Cl Der.....	PCT	999.12	4.1.12
139.04	Taurine, Pre-col FMOC Der.....	PCT		
139.05	Taurine, Pre-col AQC Der	PCT		
139.99	Taurine, Misc.	PCT		
087.00	Thiabendazole, Spectrophotometric.....	PCT	966.28	5.1.55
087.01	Thiabendazole, HPLC	PCT		
087.99	Thiabendazole, Misc.	PCT		
105.00	Thiamine, HPLC	mg/LB		
105.01	Thiamine.....	mg/LB	942.23	45.1.05
105.99	Thiamine, Misc.	mg/LB		
135.00	Threonine, Post-col Ninhydrin Der	PCT	994.12	4.1.11
135.01	Threonine, Pre-col OPA Der.....	PCT		
135.02	Threonine, Post-col OPA Der	PCT		
135.03	Threonine, Pre-col Dansyl Cl Der.....	PCT		
135.04	Threonine, Pre-col FMOC Der.....	PCT		
135.05	Threonine, Pre-col AQC Der	PCT		
135.99	Threonine, Misc.	PCT		
086.00	Tiamulin, HPLC	g/Ton		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
086.99	Tiamulin, Misc.	g/Ton		
136.00	Tryptophan, Alka-Hydrol Post-col Ninhyd.....	PCT	988.15	45.4.04
136.01	Tryptophan, Alka-Hydrol Rev Phase LC UV.....	PCT		
136.02	Tryptophan, Alka-Hydrol Post-col OPA De.....	PCT		
136.03	Tryptophan, Alka-Hydrol + IS Rev Phase LC FI..	PCT	ISO13904:2005(E)	
136.99	Tryptophan, Misc.	PCT		
088.00	Tylosin, Plate	g/Ton	962.26	5.3.24
088.01	Tylosin, Turbid	g/Ton		
088.02	Tylosin, Plate w Hyd	g/Ton		
088.03	Tylosin, HPLC	g/Ton		
088.99	Tylosin, Misc.	g/Ton		
137.00	Tyrosine, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
137.01	Tyrosine, Pre-col OPA Der.....	PCT		
137.02	Tyrosine, Post-col OPA Der	PCT		
137.03	Tyrosine, Pre-col Dansyl Cl Der	PCT		
137.04	Tyrosine, Pre-col FMOC Der	PCT		
137.05	Tyrosine, Pre-col AQC Der.....	PCT		
137.99	Tyrosine, Misc.	PCT		
000.02	Urea, Colorimetric.....	PCT	967.07	4.3.02
000.99	Urea, Misc.	PCT		
007.00	Urease Act pH Chng (BA 9-58)	PCT	AOCS	
007.99	Urease, Misc.	PCT		
138.00	Valine, Post-col Ninhydrin Der.....	PCT	994.12	4.1.11
138.01	Valine, Pre-col OPA Der.....	PCT		
138.02	Valine, Post-col OPA Der	PCT		
138.03	Valine, Pre-col Dansyl Cl Der	PCT		
138.04	Valine, Pre-col FMOC Der.....	PCT		
138.05	Valine, Pre-col AQC Der.....	PCT		
138.99	Valine, Misc.	PCT		
041.00	Vanadium, ICP	PPM		
041.99	Vanadium, Misc.	PPM		
089.00	Virginiamycin, Plate	g/Ton		
089.99	Virginiamycin, Misc.	g/Ton		
106.00	Vitamin A, Color.....	KU/LB	974.29	45.1.02
106.01	Vitamin A, UV	KU/LB		
106.02	Vitamin A, HPLC.....	KU/LB		
106.99	Vitamin A, Misc.	KU/LB		
107.00	Vitamin B12	mcg/LB	952.20	45.2.02
107.99	Vitamin B12, Misc.	mcg/LB		
108.00	Vitamin D3.....	KU/LB		
108.01	Vitamin D3, HPLC	KU/LB	982.29	45.1.22
108.02	Vitamin D3, HPLC	KU/LB		
108.99	Vitamin D3, Misc.	KU/LB		
109.00	Vitamin E	mg/KG	948.26	45.1.26
109.01	Vitamin E, Color	mg/KG	971.30	45.1.24
109.02	Vitamin E, HPLC.....	mg/KG		
109.99	Vitamin E, Misc.	mg/KG		
110.00	Xanthophyll.....	mg/LB	970.64	45.1.04
110.01	Xanthophyll.....	mg/LB		
110.99	Xanthophyll, Misc.	mg/LB		

Method Codes for AAFCO Feed Check Samples (2009)

AAFCO Code	Assay Method	Units	AOAC Reference*	AOAC Locator**
037.00	Zinc, Dithizone.....	PPM	941.03	3.3.18
037.01	Zinc, AA	PPM	968.08	4.8.02
037.02	Zinc, DCP	PPM		
037.03	Zinc, ICP, Dry Ash	PPM	968.08 (D)(a)	
037.04	Zinc, Em Spect	PPM	953.01	3.2.01
037.05	Zinc, ICP, Wet Ash	PPM	968.08 (D)(b)	
037.99	Zinc, Misc.	PPM		
090.00	Zoalene.....	PCT	961.26	5.1.57
090.01	Zoalene, HPLC.....	PCT		
090.99	Zoalene, Misc.	PCT		

*To the 18th Edition

**For the 18th Edition on-line