

Comparison of Different Methods for Determination of Amino Acids Contents in Food and Feed

Y. Zhang and J.M. Sido

(National Corn-to-Ethanol Research Center,, Edwardsville, IL)

S.N. Daniels

(Applied Biosystems, Framingham, MA)

R. Dickerson and L. Reimann

(Eurofins US in Des Moines,Des Moines, Iowa)

H. Hewitson and T.E. Wheat

(Waters Corporation, Milford, MA)

L. Novotny

(Olson Biochemistry Labs, South Dakota State University, Brookings SD)

J. Reuther, S. Ruiz and F. Sudradjat

(Eurofins Central Analytical Laboratories, Metairie, Louisiana)

K. Sjogren

(Minnesota Valley Testing Laboratories, New Ulm, MN)

Analytical Methods Used

- HPLC – ion exchange with ninhydrin post-column derivatization (3 labs)
- HPLC – ion exchange with o-phthaldialdehyde post-column derivatization (1 lab)
- HPLC – pre-column derivatization (ACCQ-Fluor Reagent Kit) with FLD (1 lab)
- UPLC – pre-column derivatization (6-aminoquinolyl-N-hydroxysuccinimidyl carbamate) (2 labs)
- LC/MS/MS – pre-column derivatization (iTRAQ™ reagent, an amine reactive compound) (2 labs).

Amino Acid Round Robin Design

- 15 amino acids:

Alanine, arginine, aspartic acid, glycine, glutamic acid, histidine, isoleucine, leucine, lysine, phenylalanine, proline, serine, threonine, tyrosine and valine

- Seven samples:

1 amino acid hydrolysate standard (NIST)

1 peanut butter (NIST)

5 DDGS hydrolysates (triplicate and 2 NIST-STD spiked hydrolysates)

- Performance Criteria:

Precision (DDGS hydrolysate triplicate)

Accuracy (amino acid hydrolysate standard from NIST, food hydrolysate from NIST, DDGS hydrolysate spiked with NIST standard)

Amino Acid Round Robin Design (AOCS)

15 amino acids:

Alanine, arginine, aspartic acid, glycine, glutamic acid, histidine, isoleucine, leucine, lysine, phenylalanine, proline, serine, threonine, tyrosine and valine

Five samples:

soybean, soybean meal, swine diet, poultry diet and BSA

Test Methods:

GC/MS/MS (1 lab)

UPLC – pre-column derivatization (6-aminoquinolyl-N-hydroxysuccinimidyl carbamate) (2 labs)

Test Scheme:

Test each sample three times a week in 5 consecutive weeks

Data Evaluation


- Variation by individual methods
- Variation by individual amino acids
- Variation by individual labs
- Variation by sample matrix

Amino Acid Round Robin Studies

Precision

Precision (RSD, %)

(Based on DDGS hydrolysate triplicate)



Precision (RSD, %)


(Based on DDGS hydrolysate triplicate)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	1.9* (0.3)	2.4 (0.7)	1.1 (0.5)	1.8 (0.7)
UPLC Pre-Column Deriv.	3.5 (1.3)	0.6 (0.1)		2.1 (2.1)
LC/MS/MS Pre-Column Deriv.	2.1 (1.4)	4.9 (1.2)		3.5 (2.0)
HPLC Pre-Column Deriv.	8			8

* Mean of 15 amino acids

Precision (RSD, %)

(Soy bean sample, tested in 5-week period, mean of 15 data points)



Precision (RSD, %)

(Soy bean sample, tested in 5-week period, mean of 15 data points)

Methods	Lab 1	Lab 2	Mean (of all labs)
GC-MS-MS	1.6* (0.5)		1.6 (0.5)
UPLC Pre-Column Deriv.	1.6 (0.9)	0.8 (0.7)	1.2 (0.5)

* Mean of 15 amino acids



Summary on Precision

- Precision is method dependent: HPLC post-column derivatization has the best performance (RSD < 2%) and HPLC pre-column derivatization has the highest RSD.
- Most methods investigated show satisfying precision for within batch and between batch runs.
- More data points needed to produce statistically sound conclusion.

Amino Acid Round Robin Study

Accuracy


Accuracy (Recovery, %)

(Based on NIST amino acid hydrolysate standard)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	98* (8)	104 (8)	105 (8)	102 (3)
UPLC Pre-Column Deriv.	96 (9)	105 (7)		100 (6)
LC/MS/MS Pre-Column Deriv.	93 (5)	107 (8)		100 (10)
HPLC Pre-Column Deriv.	105			105

* Mean of 15 amino acids

Amino Acid Contents in DDGS Hydrolysate



Total Amino Acid Contents in DDGS Hydrolysate


(ug / ml)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	247* (5)	282 (7)	256 (3)	261 (18)
UPLC Pre-Column Deriv.	216 (8)	289 (2)		253 (52)
LC/MS/MS Pre-Column Deriv.	253 (5)	250 (12)		252 (2)
HPLC Pre-Column Deriv.	266			266

* Mean of triplicate

Accuracy (Recovery, %)

(DDGS hydrolysate spiked I)



Accuracy (Recovery, %)


(DDGS hydrolysate spiked I)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	105* (5)	94 (8)	102 (7)	101 (6)
UPLC Pre-Column Deriv.	115 (4)	99 (6)		107 (11)
LC/MS/MS Pre-Column Deriv.	102 (3)	98 (13)		100 (2)
HPLC Pre-Column Deriv.	114			114

* Mean of 15 amino acids

Accuracy (Recovery, %)

(DDGS hydrolysate spiked II)



Accuracy (Recovery, %)

(DDGS hydrolysate spiked II)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	103* (5)	100 (8)	100 (8)	101 (2)
UPLC Pre-Column Deriv.	106 (5)	97 (8)		101 (6)
LC/MS/MS Pre-Column Deriv.	102 (6)	98 (17)		100 (3)
HPLC Pre-Column Deriv.	93			93

* Mean of 15 amino acids

Summary on Accuracy

(in testing hydrolysate)

- The three methods: HPLC post-column, UPLC and LC/MS/MS give similar performance, which have 100% recovery for most amino acids, with the HPLC post-column method being the most consistent.
- Except for LC/MS/MS, all the other methods have higher than 120% recovery for histidine.
- More data points needed to produce statistically sound conclusion.

Amino Acid Round Robin Studies

(Based on food and feed samples)


Record of Accuracy

(Recovery, %, based on NIST Peanut Butter)

Methods	Lab 1	Lab 2	Lab 3	Mean (of all labs)
HPLC Post-Column Deriv.	105 * (8)	114 (17)	106 (12)	108 (5)
UPLC Pre-Column Deriv.	98 (8)	120 (7)		98 (16)
LC/MS/MS Pre-Column Deriv.	101 (9)	99 (6)		100 (1)
HPLC Pre-Column Deriv.	102			102

* Mean of 15 amino acids

Amino Acid Contents in Soybean



Total Amino Acid Contents in Soy Bean (wt/wt, %, as-rec.)

Methods	Lab 1	Lab 2	Mean (of all labs)
GC-MS-MS	45.6* (0.5)		45.6
UPLC Pre-Column Deriv.	41.7 (0.5)	45.6 (0.2)	43.6 (2.7)

* Mean of 15 runs



Comparison on other Factors

Methods	Run Time (minutes)	Price (\$)
HPLC Post-Column Deriv.	60	180
UPLC Pre-Column Deriv.	10	100
LC/MS/MS Pre-Column Deriv.	24	200
HPLC Pre-Column Deriv.	40	250

Summary on Accuracy in testing food and feed

- The three methods: HPLC post-column, UPLC and LC/MS/MS give similar performance, which have 100% recovery for most amino acids, with the HPLC post-column method being the most consistent.
- More data points needed to produce statistically sound conclusion.